

# IV-Stand Normal care

## Material: stainless steel

I-N11121

### PRODUCT DESCRIPTION

Bottle Holder	Plastic, bottle holder
Load capacity per hook (kg)	2
Load capacity per system (kg)	8
Total load capacity (kg)	9.5
Material	Stainless steel
Castors	Twin castors
Base weight	non weighted
Type of height adjustment	screw height adjustment
Area of use	Normal Care
Arrangement of the hooks	90 ° to each other
Base diameter (mm)	round, 635
Base leg profile (mm)	25 x 25
Castor diameter (Ø in mm)	50
Colour of plastic parts	blue
Drip vessel	yes
Inner tube diameter (mm)	18
Length of outer tube (mm)	1,000
Maximum length/height (mm)	2,150
Minimum length/height (mm)	1,350
Number of hooks	4
Outer tube diameter (mm)	25

### PACKAGING UNIT: PAIR

Width Box1	0,19 m
Height Box1	0,10 m
Length Box1	1,19 m
Weight Box1	7,60 kg

### PACKAGING UNIT: UNIT

Width Box1	0,19 m
Height Box1	0,10 m
Length Box1	1,19 m
Weight Box1	4,10 kg



Stainless steel (1.4301)



Many different configurations



I-N11121

## IV-Stand Normal care



### Area of application

optimal hygienic characteristics (wipe disinfection)

very easy and fast assembly through a plug in and out system with only 1 screw connection

comfortable height adjustment through a screw-height adjustment with an integrated, surface-gently mold clamping element



### Technology

all plastic parts are light-resisting

castors with non-marking tread

integrated pull-out guard avoids separation of the adjustable tubes

safety bottle holder in accordance with DIN ISO 15375

a sealing cap prevents the IV-pole from hygienical sensitive liquids



### Economy

complete recyclability



Работайте с лучшими,  
всё остальное компромисс!

---

**8 (800) 775-10-98**

**medliga.ru**